


---

TRYPTIK®  
PL

---

BY SPINEART

CER  
VICAL



PLATE





# CONTENT

CONCEPT AND DESIGN	<b>PAGE 04</b>
IMPLANTS	<b>PAGE 05</b>
TECHNICAL FEATURES	<b>PAGE 06</b>
INSTRUMENT SET	<b>PAGE 07</b>
INSTRUMENTS	<b>PAGE 08</b>
SURGICAL TECHNIQUE	<b>PAGE 09</b>



# CONCEPT AND DESIGN

In December 2005, the TRYPTIK®<sub>CA</sub> cervical cage was the first ever Spineart device implanted. It was soon followed by a cervical plate, the TRYPTIK®<sub>PL</sub>, and a cervical modular cage, the TRYPTIK®<sub>MC</sub>, completing the range and creating a cervical fusion platform made of 3 devices for only 1 compact set of 7 instruments.

Like triptych, an art work composed of 3 sections, the core principle of the cervical fusion range stands on 3 solutions : a cage, a modular cage-plate, and a plate.

In each product development, Spineart is relentlessly driven by the same motto : Quality, Innovation, Simplicity.

## AT A GLANCE

### **DYNAMIC CONCEPT**

### **ANTI BACK-OUT SCREWS**

### **LOW PROFILE PLATE**

### **UNIVERSAL SET**



## INDICATIONS

The TRYPTIK® range is indicated in following pathologies between C3 to C7:

- Cervical hernia
- Cervicarthrosis
- Degenerative disc disease
- Traumatology



# IMPLANTS



TITANIUM TI 40 PLATE



TITANIUM ALLOY SCREWS

REFERENCES	
LENGHTS	4 HOLES
20mm	MOS-PL 04 20-S
25mm	MOS-PL 04 25-S
30mm	MOS-PL 04 30-S
LENGHTS	6 HOLES
35mm	MOS-PL 06 35-S
40mm	MOS-PL 06 40-S
45mm	MOS-PL 06 45-S
50mm	MOS-PL 06 50-S
55mm	MOS-PL 06 55-S
LENGHTS	8 HOLES
60mm	MOS-PL 08 60-S
65mm	MOS-PL 08 65-S
70mm	MOS-PL 08 70-S
LENGHTS	10 HOLES
80mm	MOS-PL 10 80-S
90mm	MOS-PL 10 90-S

REFERENCES	
LENGTH / DIAMETER	Ø4MM
L12mm	MOS-CS 40 12-S
L14mm	MOS-CS 40 14-S
L16mm	MOS-CS 40 16-S
L18mm	MOS-CS 40 18-S
	Ø4,5MM
L14mm	MOS-CS 45 14-S
L16mm	MOS-CS 45 16-S
L18mm	MOS-CS 45 18-S



# TECHNICAL FEATURES


## DYNAMIC CONCEPT

- The plate's oblong holes dynamically lock-up the device.




## ANTI BACK-OUT SCREWS

- The plate's anti back-out screws guarantee an optimal security and allow for 20° polyaxiality.



## LOW PROFILE PLATE

- The slimline of the plate (1.5mm) prevents from dysphagia.



## UNIVERSAL SET

- A unique set covers all of TRYPTIK®'s anterior cervical range.



# INSTRUMENT SET



#	DESCRIPTION	REFERENCE
01	SCREWDRIVER 2.5	MOS-IN 00 03-N
02	BONE AWL	MOS-IN 00 01-N
03	IMPLANT HOLDER	MOS-IN 00 02-N
04	CALIPER	ELL-IN 00 12-N
05	PLATE BENDER	MOS-IN 00 06-N
06	REVISION SCREWDRIVER	MOS-IN 00 10-N
INSTRUMENT CONTAINER		MOS-BX 10 01-N



# INSTRUMENTS

CALIPER

ELL-IN 00 12-N



BONE AWL

MOS-IN 00 01-N



IMPLANT HOLDER

MOS-IN 00 02-N



SCREWDRIVER 2.5

MOS-IN 00 03-N



PLATE BENDER

MOS-IN 00 06-N



REVISION SCREWDRIVER  
(OPTION)

MOS-IN 00 10-N






# SURGICAL TECHNIQUE

## STEP 1

### DISCECTOMY AND SELECTION OF THE PLATE SIZE



Remove the disc and insert a graft or a cage into the intervertebral disc space. According to the engraved markings on the plate bender, the size of the plate can be chosen by using the caliper.

INSTRUMENT	REFERENCE
CALIPER	ELL-IN 00 12-N

## STEP 2

### PLATE BENDING



If necessary, bend the plate with the plate bender to adapt it to the anatomy of the vertebral body.

INSTRUMENT	REFERENCE
PLATE BENDER	MOS-IN 00 06-N




# SURGICAL TECHNIQUE

## STEP 3

### PLATE POSITIONING

Screw the plate onto the implant holder and position it onto the vertebral bodies.




INSTRUMENT	REFERENCE
IMPLANT HOLDER	MOS-IN 00 02-N

## STEP 4

### PREPARATION OF THE SCREW SITES

Prepare the insertion holes of the cervical screws with the bone awl.




INSTRUMENT	REFERENCE
BONE AWL	MOS-IN 00 01-N



# SURGICAL TECHNIQUE

## STEP 5

### SCREW INSERTION



Tighten the cervical screws until the screw head is in contact with the plate, at this point the anti-back-out feature is effective.

If necessary (i.e. trauma) or desired after fusion, the screws can be removed with a removal instrument.

INSTRUMENT	REFERENCE
SCREWDRIVER 2,5	MOS-IN 00 03-N
REVISION SCREWDRIVER	MOS-IN 00 10-N

## FINAL CONSTRUCT



## NOTE

[illegible]

## NOTE

[illegible]





spineart.com

**SPINEART SA**

20 route de Pré-Bois  
1215 Geneva 15  
Switzerland

**SPINEART USA**

227 East 58<sup>th</sup> St. 2<sup>nd</sup> Floor  
New York, NY 10022  
United States of America

9200 Irvine Center Drive, Suite 150  
Irvine, CA 92618  
United States of America

Document not valid for the United States.

€ 1250

0414-VI ref. TRY-BR PL 31-E

S W I S S M A D E

